## INSTRUCTIONS FOR USE

### CAUTION – PLEASE READ CAREFULLY

Products covered in this document should only be used by dental specialists with experience in dental implantology and other specialties, such as dental diagnosis, planning, dental surgery, or prosthetic techniques. Products covered in this document designed for single use only must never be re-used. If reused, there is risk of product damage and possible prosthetic failure, which may cause health risks to the patient. Please read these instructions carefully before using the Azure Tapered Bone Level implants and components. All products covered in this document should be dry fitted before use to check that they fit correctly. The clinician is responsible for the correct use of the products, as both the planning and the procedures are under their control. Only dental specialists with appropriate experience and training should work with these products.

Please contact the manufacturer or distributor with any questions.

**U.S Federal law restricts these devices to be sold by or on the order of a licensed dentist / physician.**

An annual inspection of the prosthetic restoration and the screw by the dentist and/or the laboratory is recommended. If the screws are subject to unusual wear, the complete integrity of the prosthetic components should be checked. Failure to follow these instructions puts the patient at risk.

**During any intraoral use and handling, all Azure products must be secured to prevent them from being swallowed or aspirated by the patient.**

### DESCRIPTION

Azure Tapered Bone Level is an integrated implant system of endosseous dental implants with corresponding instrumentation and prosthetic parts.

#### Materials:

- **Tapered Bone Level Implants:** Titanium (grade 4) according to the ISO 5832-2 standard, with a surface topography obtained by a surface treatment consisting of sandblasting and double acid etching.
- **Prosthetic components:** Titanium ELI Ti-6Al-4V according to ASTM F136 and ISO 5832-3 standard.
- **Instruments and accessories:** Titanium ELI Ti-6Al-4V or stainless steel suitable for this use.

### INDICATIONS FOR USE

Azure Tapered Bone Level Implant System is indicated for the functional and esthetic oral rehabilitation of the upper or lower jaw of edentulous or partially edentulous patients.

### CONTRAINDICATIONS

Tapered Bone Level implantology products should not be placed in patients who are found to be medically unsuitable for the intended treatment.

Contra-indications include, but are not limited to: material allergies, metabolic or systemic disorders associated with wound and/or bone healing, previously irritated bones, use of pharmaceuticals that inhibit or alter natural bone remodeling, uncontrolled diabetes, alcoholism or drug addiction, blood clotting disorders, anticoagulant therapy, diseases with periodic use of high doses of steroids, metabolic bone disease, chemotherapy or radiotherapy, chronic periodontal inflammation, insufficient soft tissue coverage, any disorder that inhibits the patient's ability to maintain adequate daily oral hygiene, uncontrolled parafunctional habits, moderate or severe smoking, uncontrolled endocrine diseases, insufficient bone height and/or width, insufficient interarticular space, pregnancy, psychosis and vascular conditions.

Tapered Bone Level implantology products are not recommended for children under 18 years of age. Treatment of young adults (of 18 to 21 years of age) is not recommended until skeletal growth is complete and epiphyseal closure has occurred. The use of unloading splints is recommended in patients with parafunctional habits.

### WARNINGS AND PRECAUTIONS

- Products must be secured when handled intraorally. Mandatory precautions must be taken to avoid accidental aspiration/ingestion by the patient.
- As part of the dental implant rehabilitation and prior to surgery, the physician is solely responsible for conducting complete preoperative evaluation and planning of the entire treatment with the patient. Precise planning is necessary to avoid any accidental damage during osteotomy preparation and implant placement to vital anatomical structures, such as the mental nerve and inferior alveolar nerve, as this may result in anaesthesia, paraesthesia or dysaesthesia.
- Do not exceed recommended insertion torques.
- Implants and prosthetic components have not been evaluated for safety, compatibility, heating or migration in MRI environments.

### SPECIFIC APPLICATIONS AND INTENDED USES

#### Tapered Bone Level Implants

**Description:** Tapered Bone Level Implants can be used for immediate or delayed loading techniques following the extraction or loss of natural teeth to restore chewing function. Immediate loading is only recommended when adequate primary stability is achieved. Tapered Bone Level Implants are available in a range of diameters from 3.3 to 4.8mm.

**Intended use:** Tapered Bone Level Implants are intended for oral implantation to provide a support structure for connected prosthetic devices. The prosthetic restorations supported are single crowns, bridges and partial or full dentures. The 3.3 mm diameter implants are intended for use in reduced interdental spaces, where there is not enough alveolar bone for a larger diameter implant. Larger diameter implants are intended for both upper and lower jaw rehabilitations for functional and aesthetic rehabilitation in partially or totally edentulous patients.

**Warning and precautions:** The use of 3.3 mm diameter implants in posterior rehabilitation is not recommended. The 3.3 mm implant should only be used in the molar region if splinted. The use of electrosurgery is not recommended due to the conductivity of dental implants.

#### Procedure:

1. Bone bed preparation. If the surgical protocol used is flapless, add the soft tissue height to the depth of the reaming. Follow the guidelines for the use of drills and instruments regarding drilling speeds, intermittent drilling and adequate cooling.
2. Use the depth gauge to check the milling depth.
3. Remove the vial from the inside of the container.
4. Insert the implant insertion key previously installed in the handpiece directly into the dental implant connection.
5. Place the implant in the osteotomy previously performed and proceed with the insertion. It is recommended that the alveolus be full of blood clot at the moment of insertion to promote better integration of the implant. Implant placement with the handpiece should not exceed 15rpm.
6. **Never exceed the insertion torque of 35 Ncm for implant placement as it may damage implant connection.**  
**Caution:** The application of a higher insertion torque may damage the internal connection of the implant and cause fracture or necrosis in the bone bed. Implants can be tilted up to 45°, but any angulation greater than 30° should be splinted. Great care and attention must be taken to prevent the patient from swallowing or aspirating the components and/or tools used.
7. After implant placement, the surgeon must evaluate bone stability and quality to determine when to load the implant.  
 Depending on the selected surgical protocol, place a cover screw, healing abutment, temporary restoration or Multi-unit abutment with corresponding cover cap.

#### Healing Abutments

**Description:** Healing Abutments are used as an auxiliary transmucosal abutment. They are placed to cover the implant connection prior to prosthetic restoration to promote healing of the surrounding soft tissue after implant placement and inhibit epithelial downgrowth.

**Intended use:** Healing Abutments are intended for the formation of a soft tissue sulcus.

**Warnings and precautions:** Choose the correct height to ensure proper function and prevent transmission of masticatory forces. Prior to placement, ensure that the implant platform is free of tissue debris. Use gentle hand torque.

#### Cover Screws

**Description:** Cover Screws are used to cover the implant platform and prevent tissue overgrowth during the healing phase of the implant. The threaded portion of the Cover Screw fits inside the internal thread of the implant, while the head of the Cover Screw covers the top surface of the implant (the implant head).

**Intended use:** Cover Screws are intended to be temporarily connected to an endosseous dental implant to protect the implant connection interface during bone healing.

**Warnings and precautions:** During installation, make sure to align it to the axis of implant insertion, avoiding impacts and cross-threading which could potentially damage the thread.

### **Multi-Unit Abutments**

**Description:** Multi-Unit Abutments are used for prosthetic restorations prepared by dentist in a dental clinic. Multi-Unit Abutments are prefabricated abutments to be directly connected to the dental implant.

**Intended use:** Multi-Unit Abutments are intended to serve as a transepithelial auxiliary fixation in screw-retained prosthetic restorations.

**Warnings and precautions:** It is highly recommended to verify its position by means of periapical radiographs.

**Procedure:** In the case of straight Multi-Unit Abutments, place the appropriate abutment using the plastic holder that facilitates the insertion. Tighten the abutment using the appropriate wrench. For angled Multi-Unit Abutments, place the appropriate abutment in its desired location and position using the metal holder and then tighten the abutment screw using the implant system specific driver tip on a torque ratchet.

For torque recommendations refer to the Azure Screw Recommended Torque Document, the Azure Catalog, or Azure Tapered Bone Level Implant System Manual

### **Surgical Drills and Bone Taps**

**Description:** Surgical Instruments used as a part of dental implant treatments.

**Intended use:** Surgical Drills and Bone Taps are intended for the preparation of the cavities in the maxillary or the mandibular bone for the placement of dental implants.

**Warning and precautions:**

- Instruments must be maintained in good condition and care must be taken that instrumentation does not damage implants or other components.
- Due to the function and design of the drills, the drill tip is 0.4mm-1.1mm longer than the insertion depth of the implant.
- Use sharp, sterilised instruments. It is recommended not to exceed 20 drilling cycles.
- Maintain a low temperature during drilling. Intermittent drilling technique is recommended, with high levels of irrigation, using pre-cooled physiological saline solution or at room temperature. Drilling should be carried out with straight up and down vertical alternative movements, and constant irrigation is necessary so that the bone does not heat up and can cause bone necrosis.
- The drill must not be stopped while drilling is in progress.
- The drilling sequence described must be followed.
- Do not let residues dry on the instruments. Clean immediately after surgery.

**Surgical protocol:**

It is important to examine and evaluate each patient carefully to determine each particular situation, especially any bone or soft tissue deficits that may influence the final outcome. Depending on the density of the bone (type 1 = very hard bone, up to type 4 = very soft bone) different drilling protocols should be applied for the Tapered Bone Level Implants. This will make the adjustment of the implant in the alveolar preparation to the individual bone types more flexible depending on their quality and anatomical situation.

To create optimal healing conditions for the implant, the hard and soft tissues must be treated in a non-aggressive manner. The implant placement site must be carefully prepared, ensuring the elimination of any contaminants or sources of infection. The following measures are recommended to minimise excessive temperatures.

**Recommended drills:**

Before starting the surgery, make a base point on the bone using the pilot ball drill or Ø1.6" diameter drill.

Platform	Very soft bone (IV)	Soft bone (III)	Hard bone (II)	Very hard bone (I)
3.3	2.2 2.8* Profile Drill 3.3*	2.2 2.8 Profile Drill 3.3*	2.2 2.8 Profile Drill 3.3	2.2 2.8 Profile Drill 3.3 Bone Tap 3.3
4.1	2.2 2.8 3.5* Profile Drill 4.1*	2.2 2.8 3.5 Profile Drill 4.1*	2.2 2.8 3.5 Profile Drill 4.1	2.2 2.8 3.5 Profile Drill 4.1 Bone Tap 4.1
4.8	2.2 2.8 3.5 4.2* Profile Drill 4.8*	2.2 2.8 3.5 4.2 Profile Drill 4.8*	2.2 2.8 3.5 4.2 Profile Drill 4.8	2.2 2.8 3.5 4.2 Profile Drill 4.8 Bone Tap 4.8

\*Dense cortex status

Bone profilers (Bone Taps) are intended for use in the upper or lower jaw to prepare the osteotomy in situations involving dense bone to avoid excessive torque during implant insertion.

**Recommended speeds for drills:**

Diameter	Ø 2.2 mm	Ø 2.8 mm	Ø 3.5 mm	Ø 4.2 mm	Profile drill (3.3/4.1/4.8)	Bone Tap (3.3/4.1/4.8)
Speed	800 rpm	600 rpm	500 rpm	400 rpm	300 rpm	15 rpm

### **Implant Drivers**

**Description:** Instruments used as a part of dental implant treatments. Implant Drivers are used for the application of insertion torque to thread the implants into the bone.

**Intended use:** Intended for use to insert or remove dental implants during dental implant surgery.

**Warnings and precautions:** Never exceed 70 Ncm torque.

### **Depth Gauges**

**Description:** Instruments used as a part of dental implant treatments.

**Intended use:** Intended for use to verify the depth of an osteotomy during dental implant surgery.

### **Screwdrivers**

**Description:** Screwdrivers are reusable instruments which are used in conjunction with prosthetic components.

**Intended use:** Intended for use to tighten or loosen screws used to connect dental implant system components.

### **Ratchet and Adapters**

**Description:** Ratchet is a manual reusable wrench used to ensure that correct torque is achieved during manual tightening of dental implants and components. Adapters are used as an interface between the ratchet and the instrument used for this purpose.

**Intended use:** Intended for use to tighten dental implant and components and to ensure the correct torque is achieved during manual tightening.

Please refer to Azure™ Torque Indicating Ratchet IFU for more detailed instructions on the use, cleaning and sterilization of the ratchet.

### **Surgical Tray**

**Description:** The Surgical Tray is a reusable tray made from autoclavable material used to store and organize instruments and components during the surgical, restorative, and reprocessing procedures.

**Intended use:** Intended for use to organize, enclose, sterilize, transport, and store medical devices between surgical uses.

**Warnings and precautions:** Do not store wet instruments in the Surgical Tray.

Device	Reference	Material	Sterilization	Reuse	
TAPERED BONE LEVEL IMPLANTS	AITBL-IMPXXXX	Titanium (grade 4)	Sterilized using irradiation	Single use	
		(according to ISO 5832-2)			
HEALING ABUTMENTS	ARSBL-HAXXXXX	Titanium ELI Ti-6Al-4V	Non-sterile- to be sterilized in autoclave prior to use on the patient	Single use after sterilization	
		(according to ASTM F136 and ISO 5832-3)			
COVER SCREWS	STERILE: ASTBL-CSXXXX, ASTBL-CSXXXX	Titanium ELI Ti-6Al-4V	Sterilized using irradiation	Single use	
		(according to ASTM F136 and ISO 5832-3)			
	NON-STERILE: ASTBL-CSXXXX, ASTBL-CSXXXX	Titanium ELI Ti-6Al-4V	Non-sterile- to be sterilized in autoclave prior to use on the patient	Single use after sterilization	
		(according to ASTM F136 and ISO 5832-3)			
MULTI-UNIT ABUTMENTS	ARSBL-MNHXXXXXX, ARSBL-MXXXXXXX	Titanium ELI Ti-6Al-4V	Non-sterile- to be sterilized in autoclave prior to use on the patient	Single use after sterilization	
		(according to ASTM F136 and ISO 5832-3)			
SURGICAL DRILLS	ASTBL-NDXXXX, ASTBL-PIDXXXX, ASTBL-DXXXX, ASTBL-PRDXXXX	Stainless steel 465	Non-sterile- to be sterilized in autoclave prior to use on the patient	Reusable after sterilization- Up to 20 uses	-
		(according to ASTM F899-12B)			
BONE TAPS	ASTBL-TAPXX	Stainless steel AISI 420 C MOD	Non-sterile- to be sterilized in autoclave prior to use on the patient	Reusable after sterilization- Up to 20 uses	-
		(according to ASTM F899)			
IMPLANT DRIVERS	ASTBL-IDNCXXXX, ASTBL-IDRCXX	Stainless steel AISI 431 MOD	Non-sterile- to be sterilized in autoclave prior to use on the patient	Reusable after sterilization- Up to 20 uses	-
		(according to ASTM F899)			
DEPTH GAUGES	ASTBL-DGXX	Titanium ELI Ti-6Al-4V	Non-sterile- to be sterilized in autoclave prior to use on the patient	Reusable after sterilization- Up to 20 uses	-
		(according to ASTM F136 and ISO 5832-3)			
SCREWDRIVERS	ARA-DTXXXXX	Stainless steel AISI 420 MOD	Non-sterile- to be sterilized in autoclave prior to use on the patient	Reusable after sterilization- Up to 20 uses	-
		(according to ASTM F899)			
RATCHET & ADAPTERS	ASTBL-TIR	Stainless steel AISI 420 F MOD	Non-sterile- to be sterilized in autoclave prior to use on the patient	Reusable after sterilization- Up to 20 uses	-
		(according to ASTM F899)			
	ASTBL-RA, ASTBL-IDA	Stainless steel AISI 420 F MOD			
		(according to ASTM F899)			
SURGICAL TRAY	ASTBL-SURTRAY	PPSU, silicone	Non-sterile- to be sterilized in autoclave prior to use in a clinical setting	Reusable after sterilization	-

#### INTENDED USERS AND PATIENT GROUPS

Tapered Bone Level implantology products should only be used by dental specialists with experience in maxillary implantology and other specialties, such as dental diagnostics, planning, dental surgery, or prosthetic techniques. Their use is limited to dental laboratories and clinics.

Indicated for edentulous patients (totally or partially) who require oral rehabilitation by means of implant-supported prostheses. Rehabilitation can be unitary, multiple or overdentures, both in the upper and lower jaw. Its use is indicated in patients with complete dentofacial development. There are no differences in end users by age group, sex, ethnicity, family predisposition or genetic aspects. It is not indicated for use in patients without dental problems.

#### CLINICAL BENEFITS AND UNDESIRABLE SIDE EFFECTS

As a clinical benefit, patients can expect total or partial rehabilitation of the dentition, allowing them to regain proper masticatory function.

It is the surgeon's responsibility to provide the patient with any and all information about side effects, precautions and possible complications that may arise after implant surgery, as well as to complete the consent form.

#### SERIOUS INCIDENT NOTICE

For patients / users / third parties in the European Union with an identical regulatory regime (Regulation 2017/745/EU) if, due to the use of the product, a serious incident occurs, notify ZimVie Dental at

ZimVie Dental – Headquarters  
Monday - Friday 7:30AM - 7:00PM EST  
Address: 4555 Riverside Drive. Palm Beach Gardens, Florida 33410  
USA: 1-800-342-5454  
Canada: 1-800-363-1980

GLOBAL PHONE: 1-561-776-6700  
GLOBAL FAX: 1-561-776-1272  
Email: [emeacomplaints@zimvie.com](mailto:emeacomplaints@zimvie.com)

Manufacturer:  
Terrats Medical SL  
Calle Mogoda 75-99  
08210 - Barberá del Vallés -Barcelona (Spain)  
Tel + 34935646006  
Email: [medical.incidents@dessdental.com](mailto:medical.incidents@dessdental.com)

## STERILITY AND REUSABILITY INFORMATION

### Sterile products

Products supplied in sterile condition have been sterilized by irradiation and are intended for Single use. Sterile products must not be re-sterilized or reused under any circumstances. No responsibility will be accepted for re-sterilized devices, regardless of who performed the re-sterilization or what method was used.

**Warning:** Do not re-sterilize. Do not use the device after the expiration date stated on the label. Do not use the device if the packaging is damaged or has been previously opened.

**Caution:** Single use devices cannot be reprocessed. Reuse may cause damage or deterioration of product characteristics that may result in prosthetic solution inadequacies and/or other impairment of patient health, such as tissue infection.

### Non-sterile products

Single-use prosthetic components should be sterilized before use. The recommended method is steam autoclave treatment – Gravity displacement at **121 °C for 30 minutes, drying time 30 minutes** (according to ISO 17665-1 and ISO/TS 17665-2). Wait for the complete end of the drying cycle.

Reusable instruments and surgical trays should be cleaned and sterilized before each use. It is recommended to use the Surgical Tray for organizing instruments and components during cleaning/sterilization and during implant/prosthetic procedures.

The recommended sterilization method is the Fractioned Vacuum process (Pre-vacuum steam) for **at least 3 minutes at 132 °C (270 °F) up to 134 °C (273 °F), drying time 20 minutes**. The use of sterilization tokens is recommended, recording date and expiration date, in addition to periodic controls of the sterilization process by means of biological indicators.

Inspect reusable devices prior to each reuse and after sterilization to ensure that the integrity and performance of the product is maintained. Check the product for visible wear, deformation, or corrosion. Products showing these signs should be discarded.

Prior to use, sterilization, cleaning and disinfection are essential. Manual cleaning and disinfection tips:

- Never place instruments of different types of materials together.
- Brush and rinse with water the excess dirt and particles on the entire surface of the instruments for 25 to 35 seconds. Do not use metallic brushes to remove impurities.
- Use disposable syringes for cleaning instrument cavities.
- Immerse the instruments in a disinfectant bath suitable for dental material. When selecting detergents and disinfectants, make sure that they are products intended for this purpose and always follow the manufacturer's instructions.
- Do not use solutions containing ammonia, hydrogen peroxide or acidic substances, as these can damage the surface coating of the drills.
- Rinse thoroughly with water after application of the disinfectant.
- When cleaning the surgical tray, remove silicone holders carefully out of the cavity. Always use mild detergents even if cleaning is done using ultrasonic equipment. Dry the boxes with a soft cloth. Do not use abrasive sponges or metal parts that might damage the surface of the box.

**Caution:** Products supplied non-sterile should not be sterilized in their original packaging, sterilization pouches must be used. The Surgical Tray is not supplied in a sterile condition and is not intended in their own to maintain sterility. It should be used in conjunction with a legally marketed, validated sterilization wrap.

**Warning:** Products supplied non-sterile must be sterilized before use in the oral cavity. The use of a non-sterile device may result in tissue infection or infectious diseases.

## MAGNETIC RESONANCE IMAGING (MRI) SAFETY INFORMATION

**Warning:** The RF safety of the device has not been tested. The patient may only be imaged by landmarking at least 30 cm from the implant, or ensuring the implant is located outside of the RF coil. A patient with this device can be safely scanned on an MRI system under the following conditions:

Device name: Tapered Bone Level Implants, Healing Abutments, Cover Screws and Multi-Unit Abutments.

- Static magnetic field strength ( $B_0$ ) less than or equal to 3.0T.
- Special gradient magnetic fields of maximum 3000 gauss/cm (30 T/m).
- RF excitation: Circular Polarization (CP).
- For the body transmitting coil, a reference point must be marked at least 30 cm from the implant or ensure that the implant is outside the coil. Extremity T/R coils are allowed. Head T/R coil excluded.
- Normal mode of operation in the allowed imaging zone.
- Maximum whole body specific absorption rate of 2 W/kg (Normal mode of operation).
- Maximum head specific absorption rate not evaluated.
- Scan duration without specific limitations due to implant heating.

**Note:** removable restorations must be removed prior to scanning.

## STORAGE, HANDLING AND TRANSPORT

Sterile devices should be stored and transported dry in its original packaging at room temperature and should not be exposed to direct sunlight. Incorrect storage or transport may affect the sterile barrier of the product.

Non-sterile devices do not require special storage conditions but should be kept in a clean and dry place.

## DISPOSAL

Disposal of the devices should follow local regulations and environmental requirements, taking into account the different levels of contamination.

## SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The Summary of Safety and Clinical Performance (SSCP) for products belonging to the families of devices is available in the European Database on Medical Devices (EUDAMED). You can consult the document linked to the basic UDI-DI 84354572P01020101Z6 and UDI-DI 84354572Q0102F3 on the EUDAMED public website: <https://ec.europa.eu/tools/eudamed/#/screen/home>

**GLOSSARY OF SYMBOLS**

 Batch code	 Catalogue number	 Model number	 Medical Device
 Do not reuse	 Sterilised by radiation	 Expiry date	 Do not sterilize
 Do not use if the packaging is damaged	 Keep out of sunlight	 Non sterile	 Date of manufacture
 Medical device with CE marking by Notified Body 0051	 CE marking	 Caution	 See instructions for use Download link: <a href="http://azuredental.com/ifu">azuredental.com/ifu</a>
<b>Rx ONLY</b> Federal law in the USA restricts the sale of this device to physicians or via prescription.	 Distributor	 Magnetic resonance conditional	 <b>MANUFACTURER</b> Terrats Medical SL Calle Mogoda 75-99 08210 Barberá del Vallés, Barcelona (Spain) Tel: (+34) 34 935 646 006. Fax: (+34) 935 647 317 <a href="mailto:info@dess-abutments.com">info@dess-abutments.com</a>

For any information about Products distributed by ZimVie Dental, LLC, please contact your local Territory Manager.

  
ZimVie Dental, LLC  
4555 Riverside Drive  
Palm Beach Gardens, FL  
33410, USA  
+1-561-776-6700  
<http://www.zimvie.com/en/dental.html>

ZimVie Dental Global Headquarters  
Biomet 3i, LLC dba ZimVie Dental, LLC  
4555 Riverside Drive  
Palm Beach Gardens, FL 33410  
Phone: +1-561-776-6700  
Fax: +1-561-776-1272

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**IMPLANT CARD**

Tapered Bone Level Implants include an implantation/patient card to be filled in by the professional and given to the patient. This card allows to keep traceability of the implants used for the treatment.

Steps to follow:

1. Fill in the front side of the card with the patient's data, the date of surgery and your own data as a professional.
2. Remove the main label from the blister pack and attach it to the back of the card.

	Title	Description
	Patient identification	Indicates the patient's personal information
	Patient information website	Indicates a website where a patient can obtain additional information on the medical device
	Health care center or doctor	Indicates the Health Care Centre/doctor's address where medical information about the patient may be found
	Date	Indicates the date the information was entered or the Medical procedure took place

Azure products may not be available in your country. These instructions supersede all previous editions. For any detailed information about Azure products please contact your local distributor.

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