Azure™ Torque Indicating Ratchet

INSTRUCTIONS FOR USE

Read these instructions carefully before each use and keep them in a place where they are easily accessible to the user and relevant qualified personnel.

Carefully read the warning information indicated by this warning symbol. Improper use of the products may result in serious injury to the patient, user or third party.
To ensure proper condition and functionality for the intended use, the following instructions must be followed. Please note that improper handling of the instrument may adversely affect its service life and safety.

The medical devices are not supplied in a sterile condition and must be prepared and sterilized by the user prior to initial and subsequent use in accordance with the following instructions.

DELIVERED CONDITION, STORAGE OF BRAND-NEW PRODUCTS

As a general rule, the devices should be stored in a dry place and protected from dust, chemical vapors or components. The products are supplied without tension at approximately 10Ncm. This should be the general storage condition to ensure optimum performance and longer life.

PRODUCT DETAILS

Each individual part only belongs to the delivered instrument. Exchanging components is not permissible (even with identical instruments) and requires a new inspection to be performed by the manufacturer before the torque function can be used.

These Instructions for use are only applicable for the products specified below:

ASTBL-TIR

This product is a medical device and is only intended for use by trained dental specialists.
For reprocessing, the user must have sufficient qualifications that comply with legal provisions and training and hygiene specifications. The selection of the appropriate procedures in relation to the product is the responsibility of the user.

Intended use

This torque ratchet is designed to be used for the insertion and removal of screws or for the insertion and removal of implants within defined torque ranges for dental applications in implantology, osteosynthesis, surgery and prosthetics. The torque function can also be “locked”. In the locked position, higher torques can be used for both insertion and removal.

For instruments with ranges up to 80Ncm → using with loads above 100Ncm can damage the instrument.

Contraindication

Specific contraindications can only be considered in the context of procedures. It is therefore the user’s responsibility to select the most suitable methods and settings according to the individual anatomical characteristics of their patients. The torque ratchet should not be used in case of any intolerance or allergy to conventional surgical stainless steel.

Combination with instruments or other products

Adapters are available to allow this torque ratchet to be used with many different instruments. The user must ensure that the correct size is selected for the intended instrument connection. Due to the large number of possible combinations (including combinations with other manufacturers’ end instruments), please contact your distributor or contact our Customer Service Department.

When using adapters from other manufacturers, it is necessary to check their guidelines for compatibility with these instructions, at least with regard to the connection size to be used, the intended user and the preparation. We are not liable for any damage caused by combinations with other manufacturers’ products.

USE / HANDLING

Immediately prior to each use, the product should be inspected for signs of wear, or loss of function, or corrosion. It must also be correctly assembled. Damaged products or those with any of the above defects must be taken out of service immediately and must not be used in this condition.
Products with damaged sterile packaging (after processing by the user) must not be used and must be reprocessed according to these instructions.

Possible default settings

Prosthetic setting – torque function: the desired torque range can be continuously adjusted with the adjusting nut via the spring. The setting is visible on the scale on the scale sleeve.

Surgical setting – locked function: turn the adjusting nut to the reading mark ∞ (infinity). Do not overtighten.

(See Figure 1) Do not loosen the two screws **X** on the adjusting nut as this will cause the factory setting to be lost.


Figure 1

Tool change (adapter)

Retract the pin in the direction of the arrow (→) on both sides with thumb and forefinger and pull out or insert the tool (adapter). The adapter should “click” to ensure the correct insertion of the tool (see figure 1)

Correct handling of torque activation

- The pressure point for exact torque activation only is on the handle of the adjusting nut (see arrow in **Figure 2**).
- Activate only with one finger pressure.
- Do not grip the handle with thumb and forefinger for activation.
- When the set torque is reached, the scale sleeve is bent on the axis of the ratchet head. The activation takes place visibly and audible.

Do not press further after activating the torque. Damage to the ratchet or dental components could occur.

When the handle is released, the ratchet returns to the initial position.


Figure 2

MATERIAL

The product is made from the following materials:

- High grade stainless steel
- PEEK

(RE)PROCESSING

The torque ratchets are delivered in a state of cleanliness predefined by the manufacturer. They have therefore optimally adapted to the end-user handling described below. The working methods required for this have been validated by the manufacturer. In order to provide the user with information on a working procedure for each cleaning step, a reprocessing efficiency test has been carried out for a procedure set in an accredited laboratory. The validation of the procedure used for cleaning and sterilization is the responsibility of the user on site or the responsibility lies with a central sterile supply department (CSSD).

Any modification of the packaging or the product is also considered a modification of the validated delivery condition. In case of using several torque ratchets, the parts must not be interchanged. Each part belongs to its respective instrument.

Do not use metal brushes or cleaning sponges.

- Processing may only be carried out by persons trained for this purpose.
- The quality of the water to be used must be at least potable (see notes on the individual processing steps).
- The cleaning and disinfection agents used in the efficiency test are listed in these processing instructions. If other cleaning and/or disinfection media are used, they must be listed in the RKI or VAH and be compatible with the material. The PH value should be between 4.5 and 10.
- The processor is responsible for ensuring that the processing actually performed achieves the desired results with the instrument, materials and personnel employed in the processing facility. This usually requires routine validation and monitoring of the processing and equipment used.
- Please also observe the legal regulations in force in your country as well as the hygiene specifications of the medical practice or hospital. This applies in particular to the different specifications regarding effective deactivation of prions.

When selecting water qualities for instrument treatment, demineralized water should always be preferred so that corrosion-causing salt and silicate deposits are not added and are avoided or minimized as much as possible.

Transport / Place of use – Processing

The first steps of proper processing begin immediately after use on the patient.

Heavy accumulation of dirt, residues e.g. of filling material, disinfectant or other drugs must be removed before the instruments are deposited.

- Whenever and wherever possible, dry disposal (wetted closed system) should be preferred. Otherwise, disposal should be carried out according to the usual sanitary processes. The torque ratchet should be transported and disposed of in a closed container or in an airtight protective sleeve.
- In general, drying of any residues from the application should be avoided.
- Long waiting times until processing, e.g. at night or over the weekend, should be avoided for both types of disposals (<6 hours).

Cleaning and disinfection

Disinfectant and cleaning solutions with a pH value between 4.5 and 10 are provided for cleaning, follow the manufacturer's instructions for these products (e.g. for the intended use, dosage, action time, etc.).

In general, when storing parts for cleaning purposes, make sure that they touch or overlap as little as possible, in order to avoid wash shadows and to make the cleaning process as efficient as possible.

Basic principles

Depending on the possibilities, a mechanical process (washer-disinfectors) should preferably be used for cleaning and disinfection. Due to the considerably lower effectiveness and reproducibility, only a manual procedure (even with the use of an ultrasonic bath) should be used if a mechanical procedure is not available.

The preparation and pretreatment described below must be carried out in both cases.

Preparation for decontamination

Immediately after use (within 2 hours at the latest) the instruments should be cleaned of heavy accumulation of dirt. The torque ratchet must be separated into its individual parts before cleaning (irrespective of the selected cleaning type). This is possible without tools. To do so, only the adjusting nut must be unscrewed completely (see **Figure 3**).

The plastic washer must not be lost, as this would adversely affect the accuracy of the instrument. (The plastic washer should only be removed in case of visible soiling. If necessary, the washer can be removed. After cleaning, press the washer back in).


Figure 3

Pre-treatment process

Pre-cleaning should always be carried out regardless of the type of subsequent cleaning. Wash the products with cold running water (drinking water quality, <40 °C) until all visible soiling has been removed. Stubborn dirt should be removed with a soft brush. Cavities and lumens should be washed out with a pressure washer (or similar) intensively (>30 seconds) under cold running water (drinking water quality, <40 °C).

Mechanical process – Thermal disinfection

The general suitability of the instruments for effective mechanical cleaning and disinfection was certified by an independent test laboratory accredited and recognized by the ZLG (§15 (5) of the German Medical Products Act or MPG) using a Miele G7835 CD washer-disinfector (thermal disinfection, Miele & Cie. KG, Gütersloh) and the Neodisher® MediClean pre-cleaning and cleaning medium (Dr. Weigert GmbH & Co. KG, Hamburg). The described process was considered for this purpose.

Cleaning and disinfection unit and media

The following must be considered when selecting a cleaning and disinfection unit:

- The unit should in general have proven efficacy (e.g., **DGHM** [German Association for Hygiene and Microbiology] or **FDA** [Food and Drug Administration] approval / clearance / registration or CE marking in accordance with DIN EN ISO 15883),
- If possible, a proven program should be used for thermal disinfection (A_0 value >3000 or, in case of older devices, at least 5 minutes at 90°C / 194°F) (in case of chemical disinfection there is a danger of disinfectant residues on the instruments),
- The program used must be suitable for the instruments and contain sufficient wash cycles,
- Only demineralized water may be used for rinsing,
- The air used to dry must be filtered (oil-free, low in germs and particles),
- The cleaning and disinfection unit must be serviced and checked regularly.

It is essential to observe the material applications, concentrations, temperatures and treatment times, as well as the specifications for rinsing, indicated by the manufacturer of the cleaners and disinfectants.

Mechanical cleaning /disinfection (→RECOMMENDED)

Program parameters used for the test (Des-Var-TD program / Miele G7835 CD cleaning and disinfection unit):

- The parts should be placed in the mobile injection unit (E450/1) on a tray.
- 1 minute pre-cleaning (cold running water of drinking quality at <40 °C) → Drain water → 3 minutes pre-cleaning (cold running water of drinking quality at <40 °C) → Drain water.
- 10 min. cleaning at 55±5 °C with 0.2% alkaline cleaner (Neodisher® MediClean 0.2%) → Drain water.
- 1 minute rinsing with demineralized water at <40 °C → Water drain → 2 minutes rinsing with demineralized water at <40 °C → Water drain.
- Automatic disinfection >5 minutes at 92±2 °C with demineralized water.
- Automatic drying procedure at 90±2 °C of the cleaning and disinfection unit for at least 30 minutes (± 60±5 °C in the washing chamber).

(Re)processing process

- Insert the instruments in the cleaning and disinfection unit. Make sure that the instruments are not in contact with each other.
- Start the program.
- Remove the instruments from the cleaning and disinfection unit immediately after the end of the program and make sure that they are in sufficient dry state before packaging them.
- Check and pack the instruments without delay after removal.

Manual post-drying

If a manual post-drying is required, this should be performed with a lint-free cloth, and/or blowing of the lumen with sterile, oil-free pressurized air.

Manual process

The certificate of general suitability of the instruments for effective manual cleaning and disinfection has been issued by an independent testing laboratory, accredited and recognized by the ZLG (§15 (5) of the German Medical Devices Act or MPG) with the use of the cleaners and disinfectants mentioned below. For this purpose, the described process was considered.

Manual cleaning

1. Place the products in an alkaline cleaner (e.g. Neodisher® MediClean 0.5%) in an ultrasonic bath with an ultrasonic action time of 10 minutes. The maximum temperature of 40 °C must not be exceeded. The instructions of the cleaner manufacturer must be followed.
2. The products must be thoroughly cleaned afterwards with a soft brush. Wash the cavities and lumens, if present, with a pressure water gun (or similar) intensively (>30 seconds).
3. Rinse the products with running water (drinking water quality) to remove the cleaner (>15 seconds).

Manual disinfection

1. Immerse the products in a disinfectant listed by the RKI or VAH. Follow the instructions of the disinfectant manufacturer. It must be ensured that the disinfectant actually reaches all areas of the product (move the parts in the disinfectant bath or wash hidden surfaces with disinfectant by infection, without cannula).
2. The efficiency test for the process was performed with the disinfectant: 3% Korsorex plus (Bode Chemie, Hamburg) for 15 minutes.
3. Washing of the products (complete washing inside, outside and in the cavities) in demineralized water >60 seconds.

Manual drying

1. Manual drying with a lint-free disposable cloth. In order to avoid water residues in the cavities to a large extent, it is recommended to remove them by bowing them out with oil-free sterile compressed air. In order to avoid water residues in the cavities to a great extent, it is recommended to remove them by blowing with sterile pressurized air without oil.

Control

Thorough checks and function tests before and after use are the best possibility to detect and discard an instrument that is no longer suitable for operation. The working and operating areas (e.g. adapter housing and torque activation) as well as moving parts must be checked with particular care. Allow the parts to cool down to room temperature. Parts with damaged surfaces, flaking, dirt or discoloration and corrosion must be taken out of service. Put out of service instruments that are deformed or no longer functional or instruments that are otherwise damaged. Instruments that are still dirty must be cleaned and re-sterilized.

Maintenance

In case of using several torque ratchets, parts must not be interchanged. Each part belongs to its respective instrument.

- – Lightly lubricate the marked areas (see Figure 3) with instrument maintenance oil.

Please note that only instrument oils (white paraffinic oil without corrosion inhibitors or other additives) that are permitted for steam sterilization, taking into account the maximum sterilization temperature applied, and that have proven biocompatibility may be used and that the smallest possible.

Assemble the ratchet and perform the function test.

After assembly and before sterilization, the torque ratchet must be in a relaxed state at max. 10 Ncm.

Packaging

Sterilization of the products should be carried out in a suitable sterilization container. The manufacturer's certification was carried out using a double sterilization package (common in hospitals), i.e. suitable single sterile foil package can also be used. Rapid sterilization processes and sterilization of unpackaged instruments are not permitted.

Sterilization

The recommended sterilization method is the Fractioned Vacuum process (Pre-vacuum steam) for **at least 3 minutes at 132 °C (270 °F) up to 134 °C (273 °F), drying time 20 minutes.** The use of sterilization tokens is recommended, recording date and expiration date, in addition to periodic controls of the sterilization process by means of biological indicators. Inspect reusable devices prior to each reuse and after sterilization to ensure that the integrity and performance of the product is maintained. Check the product for visible wear, deformation, or corrosion. Products showing these signs should be discarded.

STORAGE

After sterilization, the products should be stored in the sterilization container dry and dust-free.

MATERIAL RESISTANCE

When selecting the cleaner and disinfectant, make sure that they do not include the following components:

- Organic oxides, minerals or oxidizers or stronger bases.
- Organic solutions (e.g. alcohols, ether, ketone, benzine)
- Oxidizing agents (e.g. hydrogen peroxide)
- Halogens (e.g. chlorine, iodine, bromine)
- Aromatic hydrocarbons/halogenated hydrocarbons

Acid brighteners and neutralizers may not be used.

All instruments may only be exposed to temperatures not exceeding 138 °C / 280 °F.

PRODUCT SERVICE LIFE

Once the set torque has been reached 5000 times, the service life is over. Frequent reprocessing has little effect on these instruments, provided that due care is taken, the product is not damaged and is fully functional. Wear and tear and damage from use usually determine the end of product life and depend on many factors, including the type, duration and frequency of use, as well as the handling, storage and transport of the instruments.

Damaged, blunt or dirty instruments cannot be used.

Any liability is excluded in case of non-compliance. The same applies for damage due to improper reprocessing or handling, e.g. excessive mechanical action, dropping, overloading, etc.

REPAIRS/MAINTENANCE

The factory setting has an accuracy in the range of $\pm 10\%$ of the set torque value. A repair and testing service for torque ratchets is offered. We recommend an annual check of the torque ratchet by the manufacturer. All warranties are void for instruments that have been repaired by third parties not expressly authorized by the manufacturer. Before shipment, the instrument must be visibly clean and sterilized according to these instructions (e.g. by means of a change indicator on the sealed sterilized foil). Otherwise, the manufacturer will return the untreated product to the sender without treatment.

DISPOSAL

If the instruments can no longer be repaired or processed, they should be disposed of in the usual way in the practice or clinic. Local regulations must be observed.

SERIOUS INCIDENT NOTICE

For patients / users / third parties in the European Union with an identical regulatory regime (Regulation 2017/745/EU) if, due to the use of the product, a serious incident occurs, notify ZimVie Dental at

ZimVie Dental – Headquarters

Monday - Friday 7:30AM - 7:00PM EST

Address: 4555 Riverside Drive. Palm Beach Gardens, Florida 33410

USA: 1-800-342-5454

Canada: 1-800-363-1980

GLOBAL PHONE: 1-561-776-6700

GLOBAL FAX: 1-561-776-1272

Email: emeacomplaints@zimvie.com

Manufacturer:

Terrats Medical SL

Calle Mogoda 75-99

08210 - Barberá del Vallés -Barcelona (Spain)

Tel + 34935646006

Email: medical.incidents@dessdental.com

GLOSSARY OF SYMBOLS

 Batch code	 Catalogue number	 Model number	 Medical Device
 Caution	 Non sterile	 Date of manufacture	 CE marking
 See instructions for use Download link: azuredental.com/ifu	<p style="text-align: center;">Rx ONLY</p> Federal law in the USA restricts the sale of this device to physicians or via prescription.	 Distributor	 MANUFACTURER Terrats Medical SL Calle Mogoda 75-99 08210 Barberá del Vallés, Barcelona (Spain) Tel: (+34) 34 935 646 006. Fax: (+34) 935 647 317 info@dess-abutments.com

For any information about Products distributed by ZimVie Dental, LLC, please contact your local Territory Manager.

ZimVie Dental, LLC
 4555 Riverside Drive
 Palm Beach Gardens, FL
 33410, USA
 +1-561-776-6700
<http://www.zimvie.com/en/dental.html>

ZimVie Dental Global Headquarters
 Biomet 3i, LLC dba ZimVie Dental, LLC
 4555 Riverside Drive
 Palm Beach Gardens, FL 33410
 Phone: +1-561-776-6700
 Fax: +1-561-776-1272

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Azure products may not be available in your country. These instructions supersede all previous editions. For any detailed information about Azure products please contact your local distributor.

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